

MODEL NO.
SA-06

ENGINEERING

ENGINEERING



ELASTICALLY COUPLED BEAM APPARATUS (DIGITAL)

Description:

Apparatus consists of a three parallel bar suspension system with elastic beam at their upper and lower ends. The upper ends of the two outer suspension rods are tied to a vertical wooden board while central suspension rod may be tied to the centre of another elastic beam supported at two outer ends only. Apparatus is supplied complete with a supporting stand.

Two Load Cell of 10 Kg. capacities with a digital Indicator is supplied with the apparatus. Digital Indicator is capable to transfer the data through USB to the computer are supplied with the apparatus. Educational software also supplied with the set up to make calculation and complete results/data.

EDUCATIONAL SYSTEM YOU CAN BUY FROM:

ENGINEERING MODELS & EQUIPMENT

Regd. Office:

Mahavir Jain Market
7 – Civil Lines, P.B. No. – 13
Roorkee – 247667 (U.K.)

Email: emodrke@gmail.com, enggmodels@bsnl.in

Factory & Works:

Adarsh Nagar
New Hardwar Road
Roorkee – 247667 (U.K.)
Phones: 01332-273721
Website: enggmod.com

AN ISO 9001:2000 & 14001:2004 COMPANY